Document:	Code:	Page:
Analytical Method	DER-10000-01	1 of 1

Title:		Valid from:
Analytical Method Determination of vitamin C in tablet and solution		
		Replaces:
		NEW
		Date of revision:
		SEP 2008
Prepared by:	Approved by:	Put into force by:
Dennis Eriksen	Dennis Eriksen	Dennis Eriksen

1. Purpose

The purpose of this analytical method is to describe the determination of vitamin C by redox tiltration.

2. Index

1.	Purpose	1
2.	Index	1
3.	Enclosures	1
4.	Principle	2
5.	Apparatus	2
6.	Reagents	2
7.	Method	2
8.	Calculation	3
9.	References	4

3. Enclosures

- 1. Print of method and example of titration curve, 3 pages
- 2. Dokumentation, 2 pages

Document:	Code:	Page:
Analytical Method	DER-10000-01	2 of 2

Title:		Valid from:
Analytical Method Determination of vitamin C in tablet and solution		
		Replaces:
		NEW
		Date of revision:
		SEP 2008
Prepared by:	Approved by:	Put into force by:
Dennis Eriksen	Dennis Eriksen	Dennis Eriksen

4. Principle

Vitamin C (ascorbic acid) in tablets and solutions is determined by redox titration.

5. Apparatus

Metrohm 721 Net Titrino, 728 stirrer and electrode 6.0431.100 or equivalent. The titrator is fully controlled by the PC-program TiNet.

6. Reagents

lodine 0,1 N, with known titer, e.g. Bie & Berntsen no. LAB00260 Sulphuric acid, dilute Ph.Eur. reagent no. 1086804.

7. Method

A. Tablets

Accurately weigh an amount of powdered material (= p g) containing about 150 mg. of ascorbic acid, into 250 ml. beaker. Add 20 ml. of sulphuric acid, dilute and 100 ml. of water.

Titrate with 0.1 N iodine determination the end-point by potentiometry.

For protein combi tablet with acceptance limit within 6 - 9 mg. An amount of 1.5 g. powdered material is weighed out. The method has to be changed to Abs.start at 1 ml. and Abs.stop at 8 ml. (Look at enclosure 1)

Document:	Code:	Page:
Analytical Method	DER-10000-01	3 of 3

Title:		Valid from:
Analytical Method Determination of vitamin C in tablet and solution		
		Replaces:
		NEW
		Date of revision:
		SEP 2008
Prepared by:	Approved by:	Put into force by:
Dennis Eriksen	Dennis Eriksen	Dennis Eriksen

B. Solution

Accurately weigh an amount of the sample (= p g) containing about 150 mg. of ascorbic acid into a 250 ml. beaker. Add 10 ml. of sulphuric acid, dilute and dilute with water to about 120 ml.

Titrate with 0.1 N iodine, dertermining the end-point by potentiometry.

8. Calculation

A. Tablets

 $\frac{n x f 8.806 x tm}{p} = \text{mg. of ascorbic acid/tablet}$

B. Solution

- $\frac{n x f 8.806 x d}{p} = \text{mg. of ascorbic acid/ml.}$
- n is the volume of 0.1 N iodine consumed (ml.)
- f is the factor (normality/0.1) of 0.1 N iodine
- tm: is the tablet mass (g/tablet) of the sample
- d: is the density (weight per ml.) of the sample

The TiNet can perform the calculations automatically. Enter the following into the program.

Id 1 = conversion factor, ie 8.806

Id 2 = factor of 0.1 N iodine

Id 3 = tablet mass (g/tablet) or density (g/ml.) depending on the sample type SmplSize = amount of sample (g)

Document:	Code:	Page:
Analytical Method	DER-10000-01	4 of 4

Title:		Valid from:
Analytical Method Determination of vitamin C in tablet and solution		
		Replaces:
		NEW
		Date of revision:
		SEP 2008
Prepared by:	Approved by:	Put into force by:
Dennis Eriksen	Dennis Eriksen	Dennis Eriksen

9. References

The method is not identical with any pharmacopoeia.